

# Record.

## Virtua® Medical Disc Publisher

The most advanced solution for recording patient studies and reports onto disc.

Codonics Virtua line of Medical Disc Publishers provide the fastest, easiest way to record and label multiple medical studies and reports to CD or DVD. Discs are recorded, labeled and ready to leave with the patient without workflow interruption. This innovative, DICOM-compliant network device offers several models to address every need and budget. To keep pace with multi-modality environments or high volume imaging applications, Virtua XR records up to 62 CDs or 34 DVDs an hour\* while Virtua C and Virtua E offer a full feature set in an ultra-compact and economic design.

Virtua features an advanced processor for receiving and managing studies, a robotic disc recorder and printer, and a user-friendly touch screen interface.\*\* The built-in printer produces brilliant, full-color disc labels that include patient demographics and the facility's address and logo for marketing. Customers can create their own custom labels or use Codonics disc label design service.

\*Based on a typical clinical study and network configuration.

\*\*Virtua C is accessed from any PC on your network.



Intuitive touch screen turns novices into experts in minutes.

### Features/Benefits

- Automatically records patient studies and reports without tying up your workstation or employee resources
- On-demand disc creation and labeling directly from modality workstations with complete anonymization of the patient data
- Smart multi-tasking enables reception to both drives for recording and labeling of different discs simultaneously
- Conveniently view studies using DICOM viewers, specialized viewers and custom OEM viewers\*\*\*
- DICOM Structured Reports (DSRs) are received and matched to a patient's study allowing the recording of the report and study on the same disc\*\*\*
- Scheduled Archive automatically records all studies to disc for backup/disaster recovery, recording a complete history of all archive activities on each disc\*\*\*
- The Smart Drive holds Virtua's configuration settings and is the basis for our acclaimed Sunrise Express "Swap" Warranty which ships a replacement the same day when your call is received by 3 p.m. Factory warranty up to five years.
- Full medical device compliance including Class 1 FDA and MDD CE, GMP/QSR, ISO13485: 2003, 60601-1 Safety, EMC/EMI (55011(B) and 60601-1-2) for Healthcare Facilities

This model/configuration of Virtua is not available for purchase or use in the United States.

\*\*\*Optional features. Not all features available on all models.



Codonics full color CD and DVD

Black and white text or hand-written discs



Virtua is ideal for easy, inexpensive distribution of the entire patient study and report on disc.